



DispensMate

Bottle-Top Dispenser

For safer and more efficient way to handle chemicals

Features

- Excellent chemical resistance, Inert and high strength components made of PTFE, FEP, BSG, PP
- Fully autoclavable at 121°C (only liquid accessible parts)
- Available in 4 dispensing volume ranges from 0.5mL to 50mL
- Easy to disassemble and clean for periodic maintenance
- Supplied with most popular adapters (S40, GL32, GL38, GL25, GL28)

Applications

DLAB bottle-top dispensers offer high precision and accuracy with ISO 8655 conformed calibration for use in areas such as pharmacology, biotechnology, healthcare, food and beverage, chemistry, environmental monitoring and academy&research.

- 1 Easy-to-use volume adjustment
- 2 Smooth piston design provides effortless dispensing experience
- 3 Easy-to-use closure cap protects user from exposure to reagent
- 4 Adapters to comfortably fit most laboratory reagent bottles

Specifications

| Cat. No. | Volume Range | Graduation | Systematic Error | | Random Error | |
|------------|--------------|------------|------------------|-----------|---------------|-----------|
| | mL | mL | μL | % | μL | % |
| 7032100101 | 0.5–5 | 0.1 | ± 25 | ± 0.5 | ± 10 | ± 0.2 |
| 7032100102 | 1.0–10 | 0.2 | ± 50 | ± 0.5 | ± 20 | ± 0.2 |
| 7032100103 | 2.5–25 | 0.5 | ± 125 | ± 0.5 | ± 50 | ± 0.2 |
| 7032100104 | 5.0–50 | 1.0 | ± 250 | ± 0.5 | ± 100 | ± 0.2 |



| Cat. No. | 17400131 | 17400132 | 17400133 | 17400134 | 17400135 | 17400037 | 17400021 | 17100073 | 17100553 | 17100696 | 17100767 |
|-------------|------------------------|------------------------|------------------------|------------------------|-----------------------|----------------------------|------------------|--------------------|---------------------|---|--|
| Accessories | GL32, Adapter, 45/32mm | GL38, Adapter, 45/38mm | GL25, Adapter, 32/25mm | GL28, Adapter, 32/28mm | S40, Adapter, 45/40mm | Reagent Bottle (Brown, 1L) | Filling tube (m) | Discharge tube (m) | GL45, adapter, 45mm | Discharge tube combination (adapter 25mL, 50mL) | Discharge tube combination (adapter 5mL, 10mL) |